

South African Wheelchairs – Product List

Image	Product Name	Description
	South African Wheelchairs - Fixed Entry	Frame: Chromed Steel Frame Size: 18" (46cm) only Castors: 8" PVC castor Rear Wheel: 24" Solid Spoke Wheels Armrests: Fixed armrest Footrests: Fixed Footrest Upholstery: PVC Upholstery Weight: Load Capacity 100kg Gross Weight: 19kg
	South African Wheelchairs – Fixed Standard	Frame: Chromed Steel Frame Size: 18" (46cm) only Castors: 8" PVC castor Rear Wheel: 24" Solid Mag Wheels Armrests: Fixed armrest Footrests: Fixed Footrest Upholstery: PVC/Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 19kg
	South African Wheelchairs -Fixed Advanced	Frame: Powder Coating Steel Frame Size: 18" (46cm) only Castors: 8" PVC castor Rear Wheel: 24" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 19kg

	<p>South African Wheelchairs - Entry Steel</p>	<p>Frame: Powder Coating Steel Frame Size: 18" (46cm) only Castors: 8" PVC castor Rear Wheel: 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 19kg</p>
	<p>South African Wheelchairs - Standard Steel</p>	<p>Frame: Powder Coating Steel Frame Size: 22" (56cm) only Castors: 8" Solid castor Rear Wheel: 24" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 125kg Gross Weight: 20kg</p>
	<p>South African Wheelchairs - Advanced Steel</p>	<p>Frame: Powder Coating Steel Frame Size: 14" / 18" / 20" Castors: 8" x 2" Solid castor Rear Wheel: 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 80kg - 90kg Gross Weight: 21kg</p>

	<p>South African Wheelchairs – Manual-Mate Entry</p>	<p>Frame: Mild Steel Powder Coated Size: 16" / 18" / 20" Castors: 8X1" Solid Castor Rear Wheel: 24X1" Solid Spoke Wheels Armrests: Flip Back armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 80kg - 100kg Gross Weight: 19kg</p>
	<p>South African Wheelchairs – Manual-Mate Standard</p>	<p>Frame: Steel Powder Coated Size: 18" (46cm) only Castors: 8X1" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 19kg</p>
	<p>South African Wheelchairs – Manual-Mate Advanced</p>	<p>Frame: Steel Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 120kg Gross Weight: 19kg</p>

	<p>South African Wheelchairs – Saturn Bariatric</p>	<p>Frame: Double Crossbar Reinforced Powder Coated Mild Steel Size: 22" Castors: 8X2" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 125kg Gross Weight: 21kg</p>
	<p>South African Wheelchairs – Invacare Bariatric</p>	<p>Frame: Heavy-duty Double Crossbar Steel Powder Coated Size: 20" / 22" / 24" / 26" / 28" / 30" Castors: 8X2" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 270kg Gross Weight: 31kg</p>
	<p>South African Wheelchairs – Jade Bariatric</p>	<p>Frame: Reinforced Powder Coated Mild Steel Size: 22" / 24" / 26" Castors: 8X2" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 145kg - 185kg Gross Weight: 30kg</p>

	<p>South African Wheelchairs – Big-Ease Bariatric</p>	<p>Frame: Steel with Double Cross Bar Size: 20" / 22" / 24" Castors: 8X2" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 130kg - 150kg Gross Weight: 27kg</p>
	<p>South African Wheelchairs – Steel-Star Heavy Duty</p>	<p>Frame: Steel with Double Cross Bar Size: 18" Castors: 8X2" Solid Castor Rear Wheel: 24X1" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 120kg Gross Weight: 25kg</p>
	<p>South African Wheelchairs – Steel-Star Tuffee</p>	<p>Frame: Mild Steel Powder Coated Size: 18" Castors: 8X2" Solid Castor Rear Wheel: 24X1.75" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 80kg - 100kg Gross Weight: 20kg</p>

	<p>South African Wheelchairs – Steel-Star Cruiser</p>	<p>Frame: Steel Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8X2" Solid Castor Rear Wheel: Quick Release 24" Solid Mag Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 110kg Gross Weight: 19kg</p>
	<p>South African Wheelchairs – Alu-Lite Half Back</p>	<p>Frame: Folding Aluminium Powder Coated Size: 16" / 18" / 20" Castors: 7X2" Solid Castor Rear Wheel: Quick Release 24X1" Solid Spoke Wheels Armrests: Flip Back armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 16kg</p>
	<p>South African Wheelchairs – Alu-Lite Saturn</p>	<p>Frame: Folding Aluminium Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8X1" Solid Castor Rear Wheel: Quick Release 24" Solid Mag Wheels Armrests: Flip Back armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 115kg Gross Weight: 14kg</p>

	<p>South African Wheelchairs – Alu-Lite Standard</p>	<p>Frame: Folding Aluminium Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8X2" Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 95kg Gross Weight: 15kg</p>
	<p>South African Wheelchairs – Alu-Lite Flip Back</p>	<p>Frame: Folding Aluminium Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8X2" Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 100kg Gross Weight: 11kg</p>
	<p>South African Wheelchairs – AeroLite Buddy</p>	<p>Frame: Folding Aluminium Powder Coated Size: 14" / 16" / 18" / 20" Castors: 8X2" Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Flip Back armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 14kg</p>

	<p>South African Wheelchairs – AeroLite Swift</p>	<p>Frame: Folding Aluminium Powder Coated Size: 16" / 18" Castors: 8X1" Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Flip Back armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 100kg Gross Weight: 15kg</p>
	<p>South African Wheelchairs – Invacare Action 2</p>	<p>Frame: Folding Aluminium Powder Coated Size: 18" / 20" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Weight: Load Capacity 120kg Gross Weight: 15kg</p>
	<p>South African Wheelchairs – Invacare Action 3</p>	<p>Frame: Folding Aluminium Powder Coated Size: 18" / 20" / 22" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 120kg Gross Weight: 14kg</p>

	<p>South African Wheelchairs – Invacare Action 4</p>	<p>Frame: Folding Aluminium Powder Coated Size: 18" / 20" / 22" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 160kg Gross Weight: 14kg</p>
	<p>South African Wheelchairs – LiteSpeed Style</p>	<p>Frame: Folding Aluminium Powder Coated Size: 16" / 18" / 20" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 120kg Gross Weight: 14kg</p>
	<p>South African Wheelchairs – LiteSpeed Rubix</p>	<p>Frame: Folding Aluminium Powder Coated Size: 16" / 18" / 20" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 120kg Gross Weight: 16kg</p>

	<p>South African Wheelchairs – LiteSpeed Ultra</p>	<p>Frame: Folding Aluminium Powder Coated Size: 16" / 18" / 20" Castors: 8X1Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Reinforced Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 150kg Gross Weight: 14kg</p>
	<p>South African Wheelchairs – Semi-Recliner</p>	<p>Frame: Semi-reclining Steel Powder Coated Size: 16" / 18" / 20" Castors: 8X2 Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Weight: Load Capacity 110kg Gross Weight: 24kg</p>
	<p>South African Wheelchairs – Full Recliner</p>	<p>Frame: Full Recliner Steel Powder Coated Size: 16" / 18" / 20" Castors: 8X2 Solid Castor Rear Wheel: Quick Release 24" Solid Spoke Wheels Armrests: Detachable armrest Footrests: Detachable footrest Upholstery: Nylon Upholstery Backrest: Tension Adjustable Weight: Load Capacity 100kg Gross Weight: 14kg</p>



South African Wheelchairs – Full Recliner
Power-Star

Frame: Full Recliner Steel Powder Coated
Size: 16" / 20"
Castors: 8X2" Solid Castor
Rear Wheel: Quick Release 24" Solid Mag Wheels
Armrests: Detachable armrest
Footrests: Detachable footrest
Upholstery: Vinyl Upholstery
5 Recliner Adjustments
Weight: Load Capacity 110kg
Gross Weight: 23kg